Hot Plate Equipment Description

Equipment ID:	
Manufacturer:	
Model #:	
Serial #:	
Storage Location:	
Description:	Hot Plate
	Temperature Range:
Primary Purpose:	Drying samples
Calibration Item:	Verify temperature readings
Calibration Procedure:	In-House Procedure
	Outside Calibration
	Other
Calibration Equipment:	Laboratory Surface Contact Thermometer
Calibration Interval:	12 months
Specifications:	
Specified Teleranee	As stated in the applicable specification.

