Hot Plate Calibration Worksheet

Equipment ID:				Date:		
Manufacturer:				Performed By:		
Model			Last Calibration:			
Serial #:			Nex	kt Calibration Due:		
Storage Location:						
Calibration Item:		Verify temperature readings				
Calibration Procedure:		In-House Procedure for Hot Plates				
Calibration Equipment: Calibrated surface contact thermometer readable to at least 1 °C (to at least 1 °C (2 °F)	
		Thermome	ter I.D.			
☐ Constant Temperature				Tolerance ±		
□ Temperature Range				Tolerance ±		
	I					
Target Temperature	Temperature R			js .	Hot Plate	
		1	2	3	Setting	
Comments:						
Pass / Fail						
Initial By						
					·	

