

Laboratory Thermometer Calibration Worksheet

Equipment ID: _____	Date: _____
Manufacturer: _____	Performed By: _____
Model #: _____	Last Calibration: _____
Serial #: _____	Next Calibration Due: _____
Storage Location: _____	Temperature Range: _____
Readable to: _____	

Calibration Item: Verify temperature readings at a minimum of three points
 Calibration Procedure: In-House Procedure for Verifying Laboratory Thermometers
 Calibration Equipment: NIST Traceable Ref. Thermometer
 Thermometer I.D. _____

Typical Range _____ Tolerance \pm _____

Reference Thermometer Reading	Lab Thermometer Reading	Difference	Pass / Fail
			<input type="checkbox"/> Pass <input type="checkbox"/> Fail
			<input type="checkbox"/> Pass <input type="checkbox"/> Fail
			<input type="checkbox"/> Pass <input type="checkbox"/> Fail
			<input type="checkbox"/> Pass <input type="checkbox"/> Fail
			<input type="checkbox"/> Pass <input type="checkbox"/> Fail

Pass / Fail _____
 Initial By _____