

# Oven Calibration Worksheet

Equipment ID: \_\_\_\_\_ Date: \_\_\_\_\_  
 Manufacturer: \_\_\_\_\_ Performed By: \_\_\_\_\_  
 Model #: \_\_\_\_\_ Last Calibration: \_\_\_\_\_  
 Serial #: \_\_\_\_\_ Next Calibration Due: \_\_\_\_\_  
 Storage Location: \_\_\_\_\_

Calibration Item: Verify temperature readings  
 Calibration Procedure: In-House Procedure for Verifying Ovens  
 Calibration Equipment: NIST traceable thermometers readable to at least 1 °C (2 °F)  
 Thermometer I.D. \_\_\_\_\_  
 Thermometer I.D. \_\_\_\_\_  
 Thermometer I.D. \_\_\_\_\_

Constant Temperature \_\_\_\_\_ Tolerance ± \_\_\_\_\_  
 Temperature Range \_\_\_\_\_ Tolerance ± \_\_\_\_\_

Target Temperature	Temperature Readings			Oven Setting
	1	2	3	

Comments: \_\_\_\_\_  
 \_\_\_\_\_

Pass / Fail \_\_\_\_\_  
 Initial By \_\_\_\_\_